

Nurse Call System



YC-A Single-way nurse call system



YC-B Digital indication single-way nurse call system

Calling made convenient

Now Patients Comfort.....Just A Call Away.....

Nurse Call System enables nurse to provide prompt & effective responses to patients calls at all times.

Nurse call system has adopted the most advanced microcomputer special chip, the communication is exact, prompt, unintermittent calling. It is newly developed to meet the needs of hospitals, which is one of the most advanced Nurse Call Systems in China, has the features of easy installment and operation, various and stable functions.

The Nurse call system is composed of master unit (desk top or wall on), hand set, bed unit, emergency unit, corridor display screen, door lamp.

Thanks to unique advanced microprocessor based technology because of which is designed to work for 24hrs, 365days so that you can ensure unconditional support & comfort to your patient.

Whenever patient seek nurse assistance he presses the call button provided on the handset. Particular room no. is display at the nurse counter. The indication voice sounds simultaneously. The call can also be acknowledged in the corridor display screen. This will draw the attention of nurse very soon.

YC-A Single-way nurse call system : unintermittent calling, calling storage, the indication light on the master unit will indicate the bed number with sound alarm, corridor display screen is optional.

YC-B Digital indication single-way nurse call system : on the based of YC-A, it add the digital indication function to indicate the multiple calling numbers and bed numbers, can display the time, calendar, function, self-exam, scan, door lamp etc. and also provide several indication music for your choice. Corridor display screen with clock is available.

YC-C Double-way nurse call system : On the based of YC-B, it add the functions of interphone, hands-free, audiomonitor, and can connect to the outer telephone.

YC-D Multi-function & double-way nurse call system: On the based of YC-C, it add the function of special nurse, memo, broadcast.



YC-C Double-way nurse call system




YC-D Multi-function & double-way nurse call system



I type corridor indication screen with clock




II type corridor indication screen with clock



Installed extension

Indication Picture



YC-C Double-way nurse call system