



HX-C01 家用护理床	HX-C01 Home care bed
• 规格: (mm) 1960×760×560; 床面: 1900×700 (含床母)	SPECIFICATIONS: (mm) 1960×760×560; Bed board: 1900×700(including outer bed frame)
• 功能: 三摇	FUNCTION: three shakes
• 标准配置: 4" 万向轮 (对角刹车); 输液架一支	STANDARD CONFIGURATION: one dia.4" flexible wheel(diagonal brakes); one infusion support
• 可选: 不锈钢安全护栏; 用餐板; 塑料便盆; 床垫颜色	OPTIONAL: safe guard rail with stainless steel, a board for eating, plastic bedpan, colour of mattress



HX-C02 家用护理床	HX-C02 Home care bed
• 规格: (mm) 1960×760×580; 床面: 1900×700 (含床母)	SPECIFICATIONS: (mm) 1960×760×580; Bed board: 1900×700(including outer bed frame)
• 功能: 两摇	FUNCTION: two shakes
• 标准配置: 4" 万向轮 (对角刹车); 输液架一支; 豪华欧式护栏一对	STANDARD CONFIGURATION: one dia.4" flexible wheel(diagonal brakes); one infusion support; Geminous guard rail of luxury European
• 可选: 不锈钢护栏; 床垫颜色: 灰色, 咖啡色,	OPTIONAL: stainless steel geminous guard rail; the colour of mattress: gray, brown



HX-C03 家用护理床	HX-C03 Home care bed
• 规格: (mm) 2000×900×530	SPECIFICATIONS: (mm) 2000×900×530



HX-C04 豪华双摇三折体位护理床 (带手拉便位)	HX-C04 Luxury two shakes three-folded nursing bed (with manual bedpan)
• 规格: (mm) 2120×980×540 床面: 1900×850	SPECIFICATIONS: (mm) 2120×980×540 Bed board: 1900×850
• 功能: 背部倾斜80°±5°; 腿部倾斜35°±5°;	FUNCTION: back rest 80°±5°; knee board: 35°±5°
• 材质: 床架、床面板采用冷轧钢制造, 环保抗 菌粉体涂装。ABS强化塑料床头尾板。	MATERIAL: The frame and the surface of the bed are cool roll steel, coated with anti-rust powder. ABS plastic cover base.
• 标准配置: 安全限位摇杆二支; 杂物架一个; 5寸对角刹车豪华轮; 输液杆一支; 标准铝合金护栏一对。	STANDARD CONFIGURATION: two stainless steel safety rockers; one storage basket; dia.5" luxury flexible wheel(diagonal brakes); infusion support; a couple of standard aluminum safety guards.