



HX-G01 标准对数视力灯"E"

HX-G01 The standard logarithmic visual acuity Light "E"



HX-G02 不锈钢五官科椅

·规格: (mm) 640×700×1080
·座高: 450-550

HX-G02 Stainless steel

SPECIFICATIONS: (mm) 640×700×1080
HIGH SEATS: 450-550



HX-G03 "X"光胶片观片灯

·规格:
单联、双联、三联

HX-G03 X-ray film view of taillights

SPECIFICATIONS:
single joint, double, triple



HX-G07 不锈钢贮槽启盖架

·规格: (mm) 双槽650×400×900/1400
单槽400×400×900/1400

HX-G07 Stainless steel storage tank lid shelf kai

SPECIFICATIONS: (mm)
2 slots 650×400×900/1400
1 slots 400×400×900/1400



HX-G08 不锈钢换药托脚架

HX-G08 Stainless steel replacement drug mobile rack



HX-G09 不锈钢手术圆凳

·功能: 罗杆升降
·可选: 三脚, 四脚, 五脚

HX-G09 Stainless steel operation stool

FUNCTION: a lift pole
OPTIONAL: three feet, four feet and five feet can be chosen



HX-G04 不锈钢紫外线消毒车(电脑智能控制)

·规格: (mm)
550×480×860

HX-G04 UV disinfection of stainless steel lathe

SPECIFICATIONS: (mm)
550×480×860



HX-G05 不锈钢反光灯

HX-G05 Stainless steel reflector lamp



HX-G06 单孔冷光手术灯

HX-G06 Luminescence light hole surgery



HX-G10 不锈钢泡手桶架

·规格: (mm)
单桶、双桶

HX-G10 Stainless steel stands of immersion pail for soaking hands

SPECIFICATIONS: (mm)
Single barrel, double barrel



HX-G11 不锈钢自动感应洗手台

·规格: (mm)
A: 1600×600×900+360
B: 1800×600×900+360

·可选:
电池供电感应, 插电供电感应,
感应洗手液盒。

HX-G11 Stainless steel auto-sensing hand-washing station

SPECIFICATIONS: (mm)
A: 1600×600×900+360
B: 1800×600×900+360

OPTIONAL:
battery-powered sensor, lighting power supply
sensor, sensor liquid soap box